## Dear Bonnie,

I have enclosed the form requesting information update and have highlighted the Project Options that that we would need.

I have also enclosed a sample of an "arm immobilizer" which are used after facial surgery, especially cleft lip and use one pair per child. The child is sent home with them. These restraints are very expensive and difficult to fit if even found commercially. Please let me know if it would be possible to obtain these from Care Wear volunteers.

We are so appreciative of this wonderful service that your volunteers provide. It is truly a blessing to the families receiving them.

Many thanks! Sillie Jo DRM IMMOBILIZER PATTERN (1) YORK HOSPITAL 1001 SOUTH GEORGE STREET P.O. Box 15198 YORK PA 17405-7198 COTTON BROSDCLOTH 19" · WHITE ATTN VOLUNTEERS COORDINATOR . PASTEL (BILLIE TO BARRETT) · PRINTS (717) 851-2345 WASHABLE 2) Hem each edge Turn 1/4" & turn a second time to conceal edges NON-STRETCH SEW securely by machine. 3) Fold up 7" to areate a pocket (STITCH along each side) ←<u>II'/a"</u>-(4)10%" SEW EIGHT ROWS OF STITCHES -- HOLDING THE 7" FOLDED PIECE TO (NINE SLOTS -- APPROX 1/6" width) THE BACK (18" length) \*
Cut 3 pieces of TWILL TAPE (COTTON TWILL --(5) 14"-12" WIOTH). FOLD IN half & catch in stitching @ 12", 334", & 52"

from the bottom fold.

OR STURDY RIBBON

33/4° 5么"